



Vivid S60N



Product Description

The Vivid™ S60N combines the proven breadth, quality and performance of the Vivid product line with a new and innovative software image processing platform: cSound™. The Vivid S60N is GE cardiovascular ultrasound's high-end scanner.

The cSound architecture benefits all Vivid S60N probes and applications. The Vivid S60N supports the following applications: Adult and pediatric cardiac, coronary, pediatric, neonatal, fetal heart, obstetrics, gynecology, abdominal, small parts, thyroid, adult and neonatal cephalic, peripheral vascular, musculoskeletal conventional and superficial, nerves, urology, intraoperative, intracardiac and intraluminal.

System Architecture

GE's innovative, programmable and flexible software beamforming technology, cSound, provides exceptional image quality and power compared

to conventional GE hardware-based beam forming technology. In 2D, cSound offers true confocal imaging without the limitation of focal zones or sacrifice of frame rate and spatial resolution. Using both coherent and harmonic image processing, the system provides computational power, ease of imaging, workflow flexibility and product upgradeability.

The Vivid S60N excels in the following areas:

Exceptional image quality on the Vivid S60N is created through the use of True Confocal Imaging. The technique is enabled by the cSound platform taking advantage of advanced software image reconstruction and state-of-the-art graphics computer technology. The Vivid S60N combines Ultra Definition Clarity filtering, elevation compound imaging (considering a wider slice for 2D imaging) with the 6VT-D probe, HD Imaging (balanced resolution, penetration and image uniformity), virtual convex (wider field-of-view in the far field) for the linear probes and virtual apex (larger field-of-view) for the FPA probes.

Probe Technology – The XDclear™ series of probes are designed to help deliver powerful and efficient sound waves, with high bandwidth and efficiency. XDclear probe technology provides impressive deep penetration and high sensitivity while maintaining high spatial resolution. The combination of Single Crystal, Acoustic Amplifier and Cool Stack technologies is the core technology of the XDclear series of probes.

Ease of use features make Vivid S60 an extremely productive 2D cardiovascular ultrasound system.

The combination of the touch screen control with conventional (tactile) buttons provide intuitive controls, helping the operator maintain focus on the patient and the ultrasound images during the exam. The touch screen can also provide alpha-numeric (A/N) keyboard entry or a dedicated A/N keyboard option can be purchased.

Ease of use for the operator is provided by the cSound technology delivering auto optimized excellent image quality with little manipulation along with automated tools like 2D Auto EF 2.0, AFI 2.0 Productivity Package, Cardiac Auto Doppler and Scan Assist Pro.

Ergonomic features include the "FlexFit" mechanism enabling continuous pivoting height adjustment of the control panel, allowing the user to adjust distance to the control panel while providing the adequate legroom for standing or sitting positions. In addition, the articulating monitor arm (horizontal and vertical), and lightweight transducers combine to make the Vivid S60N an extremely ergonomic-friendly cardiovascular ultrasound system.

Portability – The Vivid S60N's compact size and light weight, combined with a fold-down monitor, enables easy transportation and promotes scanning at the patient site. The battery option provides a transportation mode that keeps the system ready to scan within a few seconds of being connected to a power outlet.

The cSound platform takes GE's **Raw Data** to a new level. For image processing and reconstruction, the Vivid S60N utilizes more than 100 times the data compared to the Vivid S6.

Additionally, the Vivid S60N uses the proven Raw data format technology that allows for advanced processing on archived images by applying many of the same scan controls and **advanced quantitative tools** as are available during the original exam.

General Specifications

Dimensions and Weight

- Width: 54 cm, 21.4"
- Depth: 76 cm, 30.2"
- Height: 132 cm – 167 cm, 52.0" – 65.7"
- Minimum height with folded screen: 104 cm, 41"
- Weight: <73 kg, 161 lbs

Electrical Power

- Nominal input voltage: 100-240 VAC, 50/60 Hz
- Rated power consumption: 500 VA

Uninterruptible Power (option)

- Battery backup for standby
- In case of power failure or accidental shutdown, when power is restored within less than 20 minutes, the system automatically turns on instantly, maintaining exact system state prior to shutdown
- For longer periods of power interruptions, the system automatically saves data and changes into "Standby" state

Operating System

- Windows® 10

Console Design

- Five active probe ports
- ECG port
- Integrated HDD
- Multiple USB ports (front/back)
- Integrated DVD-R multi drive (optional)
- On-board storage for B/W thermal printer
- Integrated speakers for premium sound

- Four swivel wheels – three wheel brakes, one wheel direction lock
- Integrated cable management
- Easily accessible removable air filters for cleaning
- Front and rear handles
- Rear storage trays/baskets
- Hand rest

User Interface

Operator Keyboard

- Ergonomic FlexFit design with left/right swivel and up/down arm-mobility of keyboard and monitor permitting both physiological sitting or standing operation
- Touch keyboard with support for characters in 12 languages
- Drawer type A/N keyboard (option)
- Physical keyboard support for international characters in 7 languages (option)
- Ergonomic hard key layout
- Interactive back-lighting of application-specific push buttons – adjustable back-light intensity
- Integrated gel holders
- Easy-to-learn user interface with intelligent keyboard
- Dedicated rotary for overall gain for 2D-mode
- Dedicated gain rotary for M-mode, CFM or Doppler controlled by active mode
- Image manager on the touch screen for quick review of image clipboard contents and easy export of images and loops to remote archives or media

Touch Screen

- 12" ultra-high-resolution, wide screen format, color, multi-touch LCD screen
- Interactive user-configurable dynamic software menu
- Touch-panel control of 8 TGC sliders
- Touch-panel controls content can be set to routine or extended usage

LCD Monitor

- 21.5" wide screen High-Definition (HD) flicker-free LCD display
- 256 shades of gray and 16.7 million simultaneous colors available
- Articulated monitor arm
- LCD translation (independent of console)
 - 350 mm horizontal bidirectional
 - 150 mm vertical height adjustment
 - Swivel to any viewing direction
- Fold down and rotation lock mechanism for transportation
- Horizontal viewing angle of more than 170°
- Resolution: 1920 x 1080 pixels
- Manual backlight and digital brightness and contrast adjustment for excellent viewing in different ambient light conditions
- Tint adjustments
- Separate adjustment for external monitor brightness/contrast
- Adaptive video formats and resolution for external monitor
- Selection for screen area output to external monitor

System Overview

Probe Presets

- Cardiac
- Stress (optional)
- Abdominal
- Peripheral vascular
- Fetal heart
- Pediatrics
- Neonatal cephalic
- Adult cephalic
- Small parts
- Thyroid
- Musculoskeletal conventional
- Urology
- Transesophageal
- OB/GYN

- Intracardiac
- Intraoperative
- Coronary (part of QuickApps)
- LVO Contrast (accessed through QuickApps)
- Nerves

Operating Modes

- 2D tissue
- 2D color flow
- 2D angio flow
- Color M-mode
- Tissue velocity M-mode
- Continuous wave Doppler
- Tissue M-mode
- Pulsed wave Doppler
- Anatomical M-mode
- Curved anatomical M-mode
- Tissue velocity imaging
- Tissue tracking
- Tissue synchronization imaging (optional)
- Strain imaging (optional)
- Strain rate imaging (optional)
- Tissue velocity Doppler
- Blood flow imaging
- B-flow
- 2D stress (optional)
- AFI Automated Function Imaging (optional)
- Auto EF (optional)
- 2D virtual apex imaging
- Coded phase inversion
- Compound imaging
- Extended field-of-view (LOGIQ™ View)

Scanning Methods

- Electronic sector
- Electronic convex
- Electronic linear
- CW pencil

Transducer Types

- Sector phased array
- Convex array
- Linear array

- Single crystal matrix array
- 2D matrix array (option)

Peripheral Options

Internal peripherals

- USB B/W video printer with control from system (optional)

External peripherals

- Network printers
 - Color laser printer
 - Color video printer with control from system
- Encrypted USB memory stick
- Three-pedal configurable footswitch

External outputs

- DVI-D
- Ethernet – 10 Mbps, 100 Mbps, 1 Gbps electrically isolated
- Multiple USB 2.0 ports

Accessories (options)

- Interface cable for external ECG and external respiratory
- ECG adapter for DIN-type pediatrics electrode leads
- Cable storage box

Display Modes

- Live and stored display format: Full size and split screen, both with thumbnails, for still and cine
- Instant-review screen displays 12 simultaneous loops/images for a quick study review
- Selectable display configuration of duplex and triplex modes: Side-by-side or top-bottom during live, digital replay and clipboard image recall
- Single, dual and quad-screen view
- Simultaneous capability
 - 2D + PW/CW
 - 2D + CFM/TVI + PW
 - 2D + CFM + CW
 - 2D + CFM/Angio/TVI/SRI/TT/SI/TSI
 - 2D + M/AMM/CAMM
 - 2D + CFM/Angio/TVI/SRI/TT/SI/TSI + M/AMM/CAMM

- Real-time duplex or triplex mode
- Compound + M/CFM/PW
- 2D + color split screen (simultaneous mode)

- Selectable alternating modes
 - 2D or compound + PW
 - 2D + CW
 - 2D or compound + CFM/PW
 - 2D + CFM + CW
- Multi-image (split/quad screen)
 - Live and/or frozen
 - Independent cine playback
- Timeline display
 - Independent 2D (or compound) + PW/CW/M display
 - A choice of display formats with various sizes of 2D + PW/CW/M
- Top/bottom selectable format
- Side/side selectable format

Display Annotation

- Patient name: First, last and middle
- Patient ID
- Additional patient ID
- Age, sex and birth date
- Hospital name
- Date format: Two types selectable – MM/DD/YY, DD/MM/YY
- Time format: Two types selectable – 24 hours, 12 hours
- Gestational age from LMP/EDD/GA
- Probe name
- Map names
- Probe orientation
- Depth scale marker
- Focal zone markers
- Image depth
- Zoom depth
- B-mode
 - Gain
 - Imaging frequency
 - Frame averaging
- M-mode
 - Gain
 - Frequency
 - Time scale

- Doppler mode
 - Gain
 - Angle
 - Sample volume size and position
 - Wall filter
 - Velocity and/or frequency scale
 - Spectrum inversion
- Time scale
 - PRF
 - Doppler frequency
- Color flow Doppler mode
 - Frame rate
 - Sample volume size
 - Color scale
 - Power
 - Color baseline
 - Color threshold marker
 - Color gain
- Spectrum inversion
- Acoustic frame rate
- CINE gauge, image number/frame number
- Bodymarks: Multiple human anatomical structures
- Application/preset name
- Measurement results
- Operator message
- Displayed acoustic output
 - TIS: Thermal Index Soft Tissue
 - TIC: Thermal Index Cranial (Bone)
 - TIB: Thermal Index Bone
- MI: Mechanical index
- Power output in dB
- Biopsy guide line and zone
- Heart rate
- Trackball-driven annotation arrows
- Active mode display
- Stress protocol parameters
- Parameter annotation follow ASE standard
- Free text with word library
- Scan plane position indicator and probe temperature are displayed with all TEE probes
- Image orientation marker

General System Parameters

System Setup

- Pre-programmable M&A and annotation categories
- Different user presets per probe/application may be stored for quick access
- User programmable preset capability with administrator preset protection
- QuickApps: Factory and user programmable sub-preset feature that keeps 2D and geometry settings while allowing different color flow or contrast parameters
- System frequency: 0 – 25 MHz
- Factory default preset data, protected against modification
- User Interface languages: English, LA Spanish, French, German, Italian, Portuguese (European and Brazilian), Russian, Swedish, Norwegian, Danish, Dutch, Finnish
- User-defined annotations
- Body patterns
- Customized comment home position

Comprehensive User Manual Available on Board

Available through touch-panel utility page. User manual and service manual are included on a USB memory device with each system. A printed user manual is provided for countries where this is required.

CINE Memory/Image Memory

- 1 GB of cine memory stores up to 690 s (196,000 frames) in 2D Color mode and up to 3,100 s in PW Doppler, depending on probe and settings
- Selectable cine sequence for cine review
- Measurements/calculations and annotations on cine playback
- Scrolling timeline memory
- Dual-image cine display
- Quad-image cine display
- CINE gauge and cine image number display
- CINE review loop
- CINE review speed

Image Storage

- On-board database of patient information from past exams
- User-selectable ECG and time gated acquisition available on touch panel during live
- User-selectable prospective or retrospective capture in config
- Storage formats:
 - DICOM®-compressed or uncompressed, single/multi-frame, with/without raw data, storage via clipboard and/or seamlessly directly to destination device
 - Transfer/"Save As" JPEG, MPEG and AVI formats
- Storage devices:
 - USB memory stick: 16 GB
 - CD-RW storage: 700 MB (DVD option required)
 - DVD storage: -R (4.7 GB) (DVD option required)
 - Hard drive image storage: 0.5 TB
- Compare old images with current exam
- Reload of archived data sets
- Activation control of USB devices (for security)

Connectivity and DICOM

- Ethernet network connection
- Wireless network connection
- DICOM 3.0
- Verify (optional)
- Print (optional)
- Store (optional)
- Modality worklist (optional)
- Storage commitment (optional)
- Modality Performed Procedure Step (MPPS) (optional)
- Media exchange
- DICOM spooler (optional)
- DICOM query/retrieve (optional)
- Structured reporting – compatible with adult cardiac, pediatric cardiac, vascular and abdominal (optional)
- Media store of structured reporting

- InSite™ ExC capability for remote service/access
- Support of two patient IDs in DICOM (optional)
- Separate DICOM SR and image storage destinations (optional)
- Simultaneous transfer of DICOM to multiple destinations (optional)

Patient Archive

EchoPAC™/Patient Archive

- Integrated EchoPAC functionality adds connectivity and image analysis capability to scanner
- Data format fully compatible with offline EchoPAC review/reporting stations of same or newer vintage
- Instant access to ultrasound raw data provided by the system
- Advanced post-processing analysis
- Three user levels help organizing data security requirements
- E-signoff compatibility, with clear indications in patient management screens and report screen that a report was signed off, and by whom and at what time. The signed off report and exam cannot be changed. The “Diagnosing Physician” field is automatically assigned to the user that did the sign-off

Image and Data Management

- Exceptional workflow with instant access data management
- DICOM 3.0 support – see DICOM conformance statement for details
- Support for transfer of the proprietary raw data files within the DICOM standard – configurable per mode
- 2D, CFM or TVI data at maximum frame rate may be reviewed by scrolling or by running cine loops (can contain more than 1000 images for imaging modes)
- Image clipboard for stamp-size storage and review of stored images and loops
- Built-in patient archive with images/loops, patient information, measurements and reports

- DICOM-SR Standard structured reporting mechanism (option)
- Structured findings report tools support efficient text entries with direct editing of findings text, usability improvements, new configuration options and conclusion section
- User can enter normal values which are then compared to actual measurements
- Configurable HTML-based report function
- Report templates can be customized on board
- Reports can be printed, stored to the archive and exported in PDF, CHM (Compiled HTML) and TXT format
- ASE-based default text modules (English), user-customizable
- Internal archive data can be exported to removable image storage through DICOM media
- Internal hard disk – for storing programs, application defaults, ultrasound images and patient archive
- All data storage is based on ultrasound raw data, allowing to change gain, baseline, color maps, sweep speeds, etc., for recalled images and loops
- DICOM media – read/write images on DICOM format
- DICOM viewer embedded on media (optional and selectable in Config)
- Alphanumeric data can be exported in XML format
- JPEG export (“Save As”) for still frames
- AVI and MPEG export (“Save As”) for cineloops
- Ability to transfer Systole Only for stress echo loops to PACS

Self-contained DICOM Viewer (optional)

- Exams can be transferred to CD/DVD or USB media with an integrated “EZ DICOM CD viewer™”
- Self-contained “EZ DICOM CD viewer™” allows to review exams from media on a standard PC, without installing anything on the host PC

Tricefy® Cloud Service

- Can serve as long-term archive
- Can be used to share examinations with colleagues for information or collaboration
- Can be used to share images with patients

Insite™ Express Connection (ExC) Enables Remote Service and Training

- Easy, flexible and secure connectivity configuration. The “Contact GE” on-screen button directly generates a real-time service request to the GE online engineering or application specialist. It takes a snapshot (e.g., error logs, setup files) of the system at the time of the service request to enable analysis of problem before customer contact
- Virtual Console Observation (VCO) enables the customer to allow desktop screens to be viewed and controlled remotely over the encrypted tunnel to enable real-time training, device configuration and clinical application support
- Operation of Insite Express Connection is dependent on the infrastructure being available – check with your local GE service representative
- File transfer enables the customer (biomed or clinician) to directly transfer system information (e.g., system logs, images, parametric data) to GE product engineering teams (no patient data transferred)
- Software reload provides remote application reconstruction and recovery capabilities in the event of system corruption

Scanning Parameters

- Infinite number of effective channels
- Minimum field-of-view range (depth): 0 – 2 cm (zoom) (probe dependent)
- Maximum field-of-view range (depth): 0 – 50 cm (probe dependent)
- Width range: 10 – 120 degrees
- Continuous dynamic receive focus/continuous dynamic receive aperture
- Adjustable dynamic range, infinite upper level

- Image reverse: Right/left
- Image rotation of 0°, 180°

Tissue Imaging

General

- Variable transmit frequencies for resolution/penetration optimization
- Display zoom with zoom area control
- High-Resolution (HR) zoom – concentrates all image acquisition power into selected Region of Interest (ROI)
- Variable contour filtering – for edge enhancement
- Depth range up to 50 cm – probe specific
- Selectable grayscale parameters: Gain, reject, DDP, clarity, dynamic range and compress – can be adjusted in live, digital replay and image clipboard recall (probe dependent)
- Automatically calculated TGC curves reduces operator interaction
- Automatically calculated lateral gain
- Smart depth: automatically optimizing transmit pattern parameters according to scan-depth setting

2D Mode

- Sector tilt and width control
- Frame rate in excess of 1500 fps, depending on probe, settings and applications
- Coded octave imaging with coded phase inversion – 3rd generation harmonic tissue imaging providing improved lateral and contrast resolution over conventional fundamental imaging. Features help reduce noise, improve wall definition, and axial resolution, making it well suited for a wide variety of patient groups
- True Confocal Imaging (TCI) – ultra narrow focused two-way beam profile throughout the field-of-view, maintaining frame rate, no zone stitching, no multi-line acquisition artifacts and enhanced dynamic contrast resolution throughout field-of-view compared to conventional focal imaging

- Automatic tissue optimization – single keystroke optimizes immediately automatically and dynamically different grayscale settings with the goal of signal independent uniform gain and contrast distribution
- UD clarity and UD speckle reduce imaging – an advanced image processing technique to remove speckle in real-time examining the relative difference between neighboring pixel values and determining whether the grayscale variations have a sharp difference, follow a trend, or are random in nature
- HD imaging – real-time simultaneous acquisition at dual frequencies compounded to help reduce speckle and noise while enhancing resolution and contrast
- Multiple-angle compound imaging – multiple co-planar images from different angles combined into a single image in real-time to help enhance border definition and contrast resolution, as well as reduce angular dependence of border or edge as compared to no-compound imaging
- Elevation compounding with 6VT-D probe (built in, no user control)
- LOGIQ View: Provides the ability to construct and view a static 2D image with wider field-of-view of a given transducer. This allows viewing and measurements of anatomy that is larger than what would fit in a single image
- Virtual convex on linear probes allows a wider field-of-view in depth than without virtual convex and aims to enhance image quality
- Virtual apex provides a wider field-of-view with phased array probes, effective at certain imaging views where a wide near field may be preferred
- L/R and up/down invert, in live, digital replay or image clipboard recall
- Digital replay for retrospective review or automatic looping of images, allowing for adjustment of parameters such as gain, reject, anatomical M-mode, persistence and replay speed

- Data Dependent Processing (DDP) performs temporal processing which helps reduce random noise but leaves motion of significant tissue structures largely unaffected – can be adjusted even in digital replay
- 256 shades of gray
- Colorized 2D-mode, user-selectable in real-time, digital replay

M-mode

- Trackball steers M-mode line available with all imaging probes – max steering angle is probe dependent
- Simultaneous real-time 2D- and M-mode
- M-mode PRF 1 kHz – image data acquired is combined to give high-quality recording regardless of display scroll speed
- Digital replay for retrospective review of spectral data
- Several top-bottom formats, side-by-side format and time-motion-only format – can be adjusted in live or digital replay
- Selectable horizontal scroll speed: 1, 2, 3, 4, 6, 8, 12, 16 seconds across display
- Horizontal scroll can be adjusted in live or digital replay

Anatomical M-mode

- M-mode cursor can be adjusted at any plane
- Curved anatomical M-mode – free (curved) drawing of M-mode generated from the cursor independent from the axial plane
- Can be activated from live, digital replay or image clipboard recall
- Anatomical color and tissue velocity M-mode
- M&A capability

Color Doppler Imaging

General

- Steerable color Doppler available with all imaging probes – max steering angle is probe dependent

- Trackball-controlled ROI
- Removal of color map from the tissue during digital replay
- Digital replay for retrospective review of color or color M-mode data allowing for adjustment of parameters such as encoding principle, color priority and color gain even on stored data
- PRF settings – user-selectable
- Advanced regression wall filter gives efficient suppression of wall clutter
- For each encoding principle, multiple color maps can be selected in live and digital replay – variance maps available
- More than 65,000 simultaneous colors processed, providing a smooth display two-dimensional color maps containing a multitude of color hues
- Simultaneous display of grayscale 2D and 2D with color flow
- Color invert – user-selectable in live and digital replay
- Variable color baseline – user-selectable in live and digital replay
- Multi-variate color priority function gives delineation of disturbed flows even across bright areas of the 2D-mode image
- Color Doppler frequency can be changed independently from 2D

Color Flow Imaging

- The cSound platform with its parallel beamformer architecture allows a combination of ultra-high frame rate and increased lateral resolution compared to previous generation GE scanners
- Very high digital signal processing power, maintaining high frame rates with large ROI's even for very low PRF settings
- Frame rate in excess of 700 fps, depending on probe and settings
- Variable ROI size in width and depth
- User-selectable radial and lateral averaging to help reduce statistical uncertainty in the color velocity and variance estimates

- Data Dependent Processing (DDP) performs temporal processing and display smoothing to help reduce loss of transient events of hemodynamic significance
- Digital replay for retrospective review or automatic looping of color images, allowing for adjustment of parameters such as DDP, encoding principle, baseline shift, color maps, color priority and color gain even on frozen/recalled data
- Application-dependent, multi-variate motion discriminator helps reduce flash artifacts
- Dedicated coronary flow application
- Multiple-angle compound imaging in 2D mode is maintained while in color Doppler mode

Color Angio

- Angle-independent mode for visualization of small vessels with increased sensitivity compared to standard color flow of previous GE products

Color M-mode

- Variable ROI length and position – user-selectable
- User-selectable radial averaging to help reduce statistical uncertainty in the color velocity and variance estimates
- Selectable horizontal scroll speed: 1, 2, 3, 4, 6, 8, 12, 16 seconds across display – can be adjusted during live, digital replay or image clipboard recall
- Real-time 2D image while in color M-mode
- Same controls and functions available as in standard 2D color Doppler

Anatomical Color M-mode

- GE-patented, any plane color M-mode display derived from color Doppler cine loop
- Also applicable to tissue velocity Imaging
- M&A capability

B-flow

- B-flow is a digital imaging technique that provides real-time visualization of vascular hemodynamics by directly visualizing blood reflectors and presenting this information in a grayscale display
- Use of GE-patented techniques to boost blood echoes, and to help preferentially suppress non-moving tissue signals
- B-flow is available for most vascular and shared service applications

Blood Flow Imaging

- Combines color Doppler with grayscale speckle imaging
- Helps improve delineation of blood flow without bleeding into tissue or vessel wall

Blood Flow Angio Imaging

- Combines angio with grayscale speckle imaging

Tissue Velocity Imaging

Tissue Velocity Imaging Mode

- Myocardial Doppler imaging with color overlay on tissue image
- Tissue Doppler data can be acquired in background during regular 2D imaging
- Frame rate in excess of 1220 fps, depending on probe and settings
- The velocity of myocardial segments after entire heart cycle can be displayed in one single image
- Tissue color overlay can be removed to show just the 2D image, still retaining the tissue velocity information
- Quantitative profiles for TVI, tissue tracking, strain and strain rate can be derived
- Time markers for valve events derived from any TM mode help simplify understanding of signals in velocity traces or curved anatomical M-mode

Tissue Tracking Mode

- Real-time display of the time integral of TVI for quantitative display of myocardial systolic displacement

- Myocardial displacement is calculated and displayed as a color-coded overlay on the grayscale and M-mode image – different colors represent different displacement ranges

Tissue Synchronization Imaging Mode (option, enabled by Advanced QScan)

- Parametric imaging which gives information about synchronicity of myocardial motion
- Myocardial segments colored according to time to peak velocity, green for early and red for late peak
- Waveform trace available to obtain quantitative time to peak measurement from TSI Image
- Available in live scanning, as well as an offline calculation derived from tissue Doppler data
- Additional features in combination with multi-dimensional imaging option
- Simultaneous acquisition of tri-plane TSI images covering all standard segments in apical views
- Efficient segment specific TSI time measurements
- Immediate bulls-eye report
- Automatic calculated TSI synchrony indexes
- TSI surface mapping
- LV synchronization report template
- CRT programming protocol

Strain/Strain Rate Mode (option, enabled by Advanced QScan)

- Tissue deformation (strain) and rate of deformation (strain rate) are calculated and displayed as real-time, color-coded overlay on the 2D image
- Cine compound calculates and displays cineloops generated from a temporal averaging of multiple consecutive heart cycles
- Anatomical M-mode and curved anatomical M-mode displays (SI and SRI)

Spectral Doppler

General

- Operates in PW, HPRF and CW modes

- Trackball steerable Doppler available with all imaging probes – max steering angle is probe dependent
- Selectable Doppler frequency for enhanced optimization
- High-quality, real-time duplex or triplex operation in all Doppler modes, CW and PW, and for all velocity settings
- Frame rate control for optimized use of acquisition power between spectrum, 2D and color Doppler modes in duplex or triplex modes
- Very fast and flexible spectrum analysis with an equivalent DFT rate of 0.2 ms
- Automatic Spectrum Optimization (ASO) provides a single press, automatic, real-time optimization of PW or CW spectrum scale, and baseline display
- Dynamic gain compensation for display of flows with varying signal strengths over the cardiac cycle to help improve ease of use
- Dynamic reject gives consistent suppression of background – user-selectable in real-time, digital replay or image clipboard recall
- Digital replay for retrospective review of spectral Doppler data
- Several top-bottom formats, side-by-side format and time-motion-only format – can be adjusted in live or digital replay
- Selectable horizontal scroll speed: 1, 2, 3, 4, 6, 8, 12, 16 seconds across display – can be adjusted in live or digital replay
- Adjustable spectral Doppler display parameters: Gain, reject, compress, color maps – can be adjusted in live or digital replay
- User-adjustable baseline shift – in live, digital replay and image clipboard recall
- Adjustable velocity scale
- Wall filters with range 10-2000 Hz (velocity scale dependent)
- Angle correction with automatic adjustment of velocity scale – in live, digital replay and image clipboard recall

- Auto Doppler angle
- Stereo speakers mounted in the front panel
- Display annotations of frequency, mode, scales, Nyquist limit, wall filter setting, angle correction, acoustic power indices
- Compound in duplex

PW/HPRF Doppler

- Automatic HPRF Doppler maintains its sensitivity even for shallow depths and with the highest PRF's
- Digital velocity tracking Doppler employs processing in range and time for high-quality spectral displays
- Adjustable sample volume size of 1-16 mm (probe dependent)
- Maximum sample volume depth 30 cm

CW Doppler

- Highly sensitive steerable CW available with all phased array probes

Tissue Velocity Doppler

Physiological Traces

- Integrated three-lead ECG module
- Automatic QRS complex detection
- External ECG lead input
- Internally generated respiratory trace using ECG leads
- ECG lead selection
- Adjustable ECG QRS markers

Contrast Imaging^{1, 2}

LVO Contrast¹ (accessed through QuickApps)

- Enables contrast applications intended for imaging of the left ventricle

¹ Schering developed harmonic imaging for supporting contrast agent imaging.

² GE Healthcare's Vivid scanner is designed for compatibility with commercially available contrast agents. Because the availability of these agents is subject to government regulation and approval, product features intended for use with these agents may not be commercially marketed nor made available before the contrast agent is approved for use. Advanced contrast features are only enabled on systems for delivery in countries or regions where the agents are approved for use or for investigational or research use.

- LV contrast (3Sc-RS, 6VT-D and 6Tc-RS probes) enhances delineation of the LV border in combination with ultrasound contrast agents. The implementation of GE's Coded Phase Inversion (CPI) provides high-resolution detection of contrast in the LV cavity and excellent suppression of myocardial tissue signals
- LVO stress (3Sc-RS probe) provides enhanced delineation of the LV border when contrast is used as part of an exercise stress exam, preserving an adequately long continuous capture buffer length
- Comprehensive set of adult and pediatric cardiac measurements and calculations to help assess dimensions, flow properties and other functional parameters of the heart
- Comprehensive set of shared service measurements and calculations covering vascular, abdominal, obstetrics and other application areas
- Configuration package to set up a customized set and sequence of measurements to use, defining user-defined measurements and changing settings for the factory-defined measurements
- Stress echo support allowing wall motion scoring and automatic stress level labeling of measurements
- Support for measuring on DICOM images
- Cardiac Auto Doppler automatically provides Doppler measurement results for the most common parameters, with minimal user guidance
- Automatic Doppler trace functionality for use in non-cardiac applications in both live and replay
- Worksheet for review, edit and deletion of performed measurements
- Reporting support allowing a configurable set of measurements to be shown in the exam report
- DICOM SR export of measurement data
- On-board IMT package facilitates non-interrupted workflow – fully integrated with M&A, worksheet, archiving and reporting functions
- Algorithm provides robust, quick, reliable measurements which can be stored to the on-board archive for review and reporting
- IMT measurement can be made from frozen images or images retrieved from archive
- IMT package supports measurements of different regions of the intima in the carotid vessel (e.g., Lt./Rt./CCA/ICA etc.)
- Frame for IMT measurement can be selected in relation to the ECG waveform

Automatic Optimization

- Dynamic optimization of B-mode image to improve contrast resolution, TGC and grayscale (soft or sharp, user-selectable)
- Auto Spectral Optimize (ASO) – dynamic adjustments of baseline, and PRF (on live image) and angle correction

Measurement and Analysis (M&A)

- Personalized measurement protocols allow individual set and order of M&A items
- Measurements can be labeled seamlessly by using protocols or post assignments
- Measurements assignable to protocol capability
- Parameter annotation follow ASE standard
- Seamless data storage and report creation
- User-assignable parameters

Intima Media Thickness (IMT) Measurements (optional)

- Automatic measurements (patent pending) of carotid artery Intima-Media Thickness (IMT) on any acquired frame

Z-Scores

- Support for three sets of user-selectable Z score publications³ covering the most common pediatric dimension measurements

Quantitative Analysis Package (Q-Analysis) (optional)

- Traces for tissue velocity or derived parameters (strain rate, strain, displacement) inside defined regions of interest as function of time
- Contrast analysis with traces for grayscale intensity or angio power inside defined regions of interest as function of time
- Curved anatomical M-mode display allowing an M-mode along an arbitrary curve in a 2D image
- Sample-area points may be dynamically anchored to move with the tissue when running the cineloop
- Cine compound displays cineloops generated from a temporal averaging of multiple consecutive heart cycles

Automated Function Imaging (AFI 2.0) (optional)

- Second generation parametric imaging tool which gives quantitative data for global and segmental wall strain
- Allows comprehensive assessment at a glance by combining three longitudinal views into one comprehensive bulls-eye view

³ C Kampmann, C M Wiethoff, A Wenzel, et. al. Normal Values of M Mode Echocardiographic Measurements of More Than 2000 Healthy Infants and Children in Central Europe. *Heart* 2000; 83; 667-672.

M Cantinotti, MD; M Scalese, MS; B Murzi, MD; et. al. Echocardiographic Nomograms for Chamber Diameters and Areas in Caucasian Children. *Journal of American Society of Echocardiography* December 2014; Volume 27, Issue 12; 1279-1292.e2.

M Cantinotti, MD; M Scalese, MS; B Murzi, MD; et. al. Echocardiographic Nomograms for Ventricular, Valvular and Arterial Dimensions in Caucasian Children with a Special Focus on Neonates, Infants and Toddlers. *Journal of American Society of Echocardiography* February 2014; Volume 27, Issue 2; 179-191.e2.

- Integrated into M&A package with specialized report templates
- 2D strain based data moves into clinical practice
- Simplified and flexible workflow with fully automated ROI tracing (if configured), adaptive ROI width and combined display of traces from all segments
- User-selectable endo or full wall global strain values displayed
- Random sequence of analysis of the three views supported
- Ability to exit tool after one or two views completed ("Easy AFI," only global strain supported)
- Applicable to transthoracic 2D and to TEE data
- Integrated AutoEF calculation

Automated Ejection-Fraction Calculation (AutoEF 2.0) (optional)

- Second generation automated EF measurement tool based on 2D-speckle tracking algorithm and on Simpson
- Performed on apical 4 Chamber and/or apical 2 Chamber views, in any order
- Integrated into M&A package with worksheet summary

Generic Measurements

- BSA (Body Surface Area)
- MaxPG (Maximum Pressure Gradient)
- MeanPG (Mean Pressure Gradient)
- % Stenosis (Stenosis Ratio)
- PI (Pulsatility Index)
- RI (Resistivity Index)
- HR (Heart Rate) – beats/minute
- A/B Ratio (Velocities Ratio)
- TAMAX (Time Averaged Maximum Velocity) – Trace method is Peak or Manual
- TAMIN (Time Averaged Minimum Velocity) – Trace method is Floor
- TAMEAN (Time Averaged Mean Velocity) – Trace method is Mean
- Volume

OB/GYN Application Module

- OB package for fetal growth analysis containing more than 100 biometry tables
- Dedicated OB/GYN reports
- Fetal graphical growth charts
- Growth percentiles
- Multi-gestational calculations (up to four)
- Programmable OB tables
- Expanded worksheets
- User-selectable fetal growth parameters based on European, American or Asian methods charts
- GYN package for ovary and uterus measurements and reporting

OB Measurements/Calculations

- Gestational age by:
 - GS (Gestational Sac)
 - CRL (Crown Rump Length)
 - FL (Femur Length)
 - BPD (Biparietal Diameter)
 - AC (Abdominal Circumference)
 - HC (Head Circumference)
 - APTD x TTD (Anterior/Posterior Trunk Diameter by Transverse Trunk Diameter)
 - LV (Length of Vertebra)
 - FTA (Fetal Trunk Cross-sectional Area)
 - HL (Humerus Length)
 - BD (Binocular Distance)
 - FT (Foot Length)
 - OFD (Occipital Frontal Diameter)
 - TAD (Transverse Abdominal Diameter)
 - TCD (Transverse Cerebellum Diameter)
 - THD (Thorax Transverse Diameter)
 - TIB (Tibia Length)
 - ULNA (Ulna Length)
- Estimated Fetal Weight (EFW) by:
 - AC, BPD
 - AC, BPD, FL
 - AC, BPD, FL, HC
 - AC, FL
 - AC, FL, HC
 - AC, HC
 - EFBW

- Calculations and Ratios
 - FL/BPD
 - FL/AC
 - FL/HC
 - HC/AC
 - CI (Cephalic Index)
 - AFI (Amniotic Fluid Index)
 - CTAR (Cardio-Thoracic Area Ratio)
- Measurements/calculations by: ASUM, ASUM 2001, Berkowitz, Bertagnoli, Brenner, Campbell, CFEF, Chitty, Eik-Nes, Ericksen, Goldstein, Hadlock, Hansmann, Hellman, Hill, Hohler, Jeanty, JSUM, Kurtz, Mayden, Mercer, Merz, Moore, Nelson, Osaka University, Paris, Rempen, Robinson, Shepard, Shepard/Warsoff, Tokyo University, Tokyo/Shinozuka, Yarkoni
- Fetal graphical trending
- Growth percentiles
- Multi-gestational calculations (four)
- Fetal qualitative description (anatomical survey)
- Fetal environmental description (biophysical profile)
- Programmable OB tables
- Over 20 selectable OB calculations
- Expanded worksheets

GYN Measurements/Calculations

- Right ovary length, width, height
- Left ovary length, width, height
- Uterus length, width, height
- Cervix length, trace
- Ovarian volume
- ENDO (endometrial thickness)
- Ovarian RI
- Uterine RI
- Follicular measurements
- Summary reports

Abdominal Calculations

- Splenic index
- Liver volume, mass, cyst
- Pancreas
- CBD
- GB wall, length

- Aorta prox, mid, dist
- Aorta iliac
- Spleen volume
- Bladder, post void bladder volume
- Renal
- Cortex thickness
- Mesenteric (CA, SMA, IMA)

Vascular Calculations

- RT ECA (Right External Carotid Artery Velocity)
- RT CCA (Right Common Carotid Artery Velocity)
- RT BIFURC (Right Carotid Bifurcation Velocity)
- RT ICA (Right Internal Carotid Artery Velocity)
- RT ICA/CCA (Right Internal Carotid Artery Velocity/Common Carotid Artery Velocity Ratio)
- LT ECA, LT CCA, LT BIFURC, LT ICA, LT ICA/CCA (same as above, for Left Carotid Artery)
- RT BULB (Right Bulbus Artery), RT VERT (Right Vertebral Artery), RT SUBC (Right Subclavian Artery), RT INN (Right Inn Artery)
- LT BULB, LT VERT, LT SUBC, LT INN
- Stent, pre-stent, post-stent
- A/B Ratio (Velocities Ratio)
- % Stenosis (Stenosis Ratio)
- S/D Ratio (Systolic Velocity/Diastolic Velocities Ratio)
- PI (Pulsatility Index)
- RI (Resistivity Index)
- HR (Heart Rate) – beats/minute
- UEV (Upper Extremity Vein velocities): IJV, SUBC, Axill V, Bas V, RV, UV, Ves, Pseudo, AVF, CephV
- UEA (Upper Extremity Artery velocities): Inn, SUBC, Axill, BA, RA, UA, Pseudo, AVF, Ves
- LEV (Lower Extremity Vein velocities): CFV, Saph FemJunc V, PopV, PTV, ATV, FV, GSV Calf, GSV Thigh, GSV Access, LSV, Saph PopJunc
- LEA (Lower Extremity Artery velocities): EIA, SFA, Pop, PTA, Peron, DPA, ATA, CFA, DFALEA

- MCA (Middle Cerebral Artery), ACA (Anterior Cerebral Artery), PCA (Posterior Cerebral Artery), AcomA (Anterior Communicating Artery), PComA (Posterior Communicating Artery), Basilar (Basilar Artery), Ves

Cardiac Measurements

- %FS (LV Fractional Shortening)
- %IVS Thck (IVS Fractional Shortening)
- %LVPW Thck (LV Posterior Wall Fractional Shortening)
- Ao Arch Diam (Aortic Arch Diameter)
- Ao asc (Ascending Aortic Diameter)
- Ao Desc Diam (Descending Aortic Diameter)
- Ao Isthmus (Aortic Isthmus)
- Ao Root Diam (Aortic Root Diameter)
- AR ERO (PISA: Regurgitant Orifice Area)
- AR Flow (PISA: Regurgitant Flow)
- AR PHT (AV Insuf. Pressure Half Time)
- AR Rad (PISA: Radius of Aliased Point)
- AR RF (Regurgitant Fraction over the Aortic Valve)
- AR RV (PISA: Regurgitant Volume Flow)
- AR Vel (PISA: Aliased Velocity)
- AR Vmax (Aortic Insuf. Peak Velocity)
- AR VTI (Aortic Insuf. Velocity Time Integral)
- ARed max PG (Aortic Insuf. End-Diastole Pressure Gradient)
- ARed Vmax (Aortic Insuf. End-Diastolic Velocity)
- AV Acc Slope (Aortic Valve Flow Acceleration)
- AV Acc Time (Aortic Valve Acceleration Time)
- AV AccT/ET (AV Acceleration to Ejection Time Ratio)
- AV EOA I (VTI) (Aortic Valve Effective Orifice Area Index by Continuity Equation VTI)
- AV EOA I Vmax (Aortic Valve Effective Orifice Area Index by Continuity Equation Peak V)

- AV CO (Cardiac Output by Aortic Flow)
- AV Cusp (Aortic Valve Cusp Separation, 2D)
- AV Dec Time (Aortic Valve Deceleration Time)
- AV Diam (Aortic Diameter, 2D)
- AV max PG (Aortic Valve Peak Pressure Gradient)
- AV mean PG (Aortic Valve Mean Pressure Gradient)
- AV SV (Stroke Volume by Aortic Flow)
- AV Vmax (Aortic Valve Peak Velocity)
- AV Vmean (AV Mean Velocity)
- AV VTI (Aortic Valve Velocity Time Integral)
- AVA (Vmax) (AV Area by Continuity Equation by Peak V)
- AVA (VTI) (AV Area by Continuity Equation VTI)
- AVA Planimetry (Aortic Valve Area)
- AVET (Aortic Valve Ejection Time)
- CO (Teich) (Cardiac Output, M-mode, Teicholtz)
- D-E Excursion (MV Anterior Leaflet Excursion)
- E' Avg (Averaged early diastolic mitral valve annular velocity)
- E' Lat (Early diastolic mitral valve lateral annular velocity)
- E' Sept (Early diastolic mitral valve septal annular velocity)
- E/E' Avg (Mitral inflow E velocity to E' Avg ratio)
- E/E' Lat (Mitral inflow E velocity to E' Lat ratio)
- E/E' Sept (Mitral inflow E velocity to E' Sept ratio)
- EDV (Cube) (Left Ventricle Volume, Diastolic, 2D, Cubic)
- EF (A-L A2C) (Ejection Fraction 2CH, Single Plane, Area-Length)
- E-F Slope (Mitral Valve E-F Slope)
- EPSS (E-Point-to-Septum Separation, M-mode)
- ERO (Effective Regurgitant Orifice)
- ESV (Cube) (Left Ventricle Volume, Systolic, 2D, Cubic)

- HR (Heart Rate, 2D, Teicholtz)
- IVC (Inferior Vena Cava)
- IVCT (Isovolumic Contraction Time)
- IVRT (Isovolumic Relaxation Time)
- IVSd (Interventricular Septum Thickness, Diastolic, 2D)
- VSs (Interventricular Septum Thickness, Systolic, 2D)
- LA Diam (Left Atrium Diameter, 2D)
- LA Major (Left Atrium Major)
- LA Minor (Left Atrium Minor)
- LA/Ao (LA Diameter to AoRoot Diameter Ratio, 2D)
- LAAd (A2C) (Left Atrium Area, Apical 2C)
- LAEDV (A-L) (LA End Diastolic Volume, Area-Length)
- LAEDV Index (A-L) (LA End Diastolic Volume Index, Area-Length)
- LAESV (A-L) (LA End Systolic Volume, Area-Length)
- LAESV Index (A-L) (LA End Systolic Volume Index, Area-Length)
- LAEDV MOD (LA End Diastolic Volume MOD)
- LAESV MOD (LA End Systolic Volume MOD)
- LIMP (Left Index of Myocardial Performance)
- LVA (s) (Left Ventricular Area, Systolic, 2CH)
- LVAd (A2C) (Left Ventricular Area, Diastolic, 2CH)
- LVAd (sax) (LV area, SAX, Diastolic)
- LVAend (d) (LV Endocardial Area, SAX)
- LVAepi (d) (LV Epicardial Area, SAX)
- LVAs (A4C) (Left Ventricular Area, Systolic, 4CH)
- LVAs (sax) (LV area, SAX, Systolic)
- LVd Mass (LV Mass, Diastolic, 2D)
- LVd Mass (LV Mass, Diastolic, M-mode)
- LVd Mass Index (LV Mass Index, Diastolic, 2D)
- LVEDV (A-L A2C) (LV Volume, Diastolic, 2CH, Area-Length)
- LVESV (A-L A2C) (LV Volume, Systolic, 2CH, Area-Length)
- LVET (Left Ventricle Ejection Time)
- LVIDd (LV Internal Dimension, Diastolic, 2D)
- LVIDs (LV Internal Dimension, Systolic, 2D)
- LVLd (apical) (Left Ventricular Length, Diastolic, 2D)
- LVLs (apical) (Left Ventricular Length, Systolic, 2D)
- LVOT Area (Left Ventricle Outflow Tract Area)
- LVOT CO (Cardiac Output by Aortic Flow)
- LVOT Diam (Left Ventricular Outflow Tract Diameter)
- LVOT max PG (LVOT Peak Pressure Gradient)
- LVOT mean PG (LVOT Mean Pressure Gradient)
- LVOT SI (Stroke Volume Index by Aortic Flow)
- LVOT SV (Stroke Volume by Aortic Flow)
- LVOT Vmax (LVOT Peak Velocity)
- LVOT Vmean (LVOT Mean Velocity)
- LVOT VTI (LVOT Velocity Time Integral)
- LVPWd (Left Ventricular Posterior Wall Thickness, Diastolic, 2D)
- LVPWs (Left Ventricular Posterior Wall Thickness, Systolic, 2D)
- LVs Mass (LV Mass, Systolic, 2D)
- LVs Mass Index (LV Mass Index, Systolic, 2D)
- LAAd (A2C) (Left Atrium Area, Apical 2C)
- MCO (Mitral Valve closure to Opening)
- MP Area (Mitral Valve Prosthesis)
- MR Acc Time (MV Regurg. Flow Acceleration)
- MR ERO (PISA: Regurgitant Orifice Area)
- MR Flow (PISA: Regurgitant Flow)
- MR max PG (Mitral Regurg. Peak Pressure Gradient)
- MR Rad (PISA: Radius of Aliased Point)
- MR RF (Regurgitant fraction over the Mitral Valve)
- MR RV (PISA: Regurgitant Volume Flow)
- MR Vel (PISA: Aliased Velocity)
- MR Vmax (Mitral Regurg. Peak Velocity)
- MR Vmean (Mitral Regurg. Mean Velocity)
- MR VTI (Mitral Regurg. Velocity Time Integral)
- MV A Dur (Mitral Valve A-Wave Duration)
- MV A Velocity (MV Velocity Peak A)
- MV Acc Slope (Mitral Valve Flow Acceleration)
- MV Acc Time (Mitral Valve Acceleration Time)
- MV Acc/Dec Time (MV: Acc.Time/Decel.Time Ratio)
- MV an diam (Mitral Valve Annulus Diameter, 2D)
- MV CO (Cardiac Output by Mitral Flow)
- MV Dec Slope (Mitral Valve Flow Deceleration)
- MV Dec Time (Mitral Valve Deceleration Time)
- MV E Velocity (MV Velocity Peak E)
- MV E/A Ratio (Mitral Valve E-Peak to A-Peak Ratio)
- MV max PG (Mitral Valve Peak Pressure Gradient)
- MV mean PG (Mitral Valve Mean Pressure Gradient)
- MV PHT (Mitral Valve Pressure Half Time)
- MV Reg Frac (Mitral Valve Regurgitant Fraction)
- MV SI (Stroke Volume Index by Mitral Flow)
- MV SV (Stroke Volume by Mitral Flow)
- MV Time to Peak (Mitral Valve Time to Peak)
- MV Vmax (Mitral Valve Peak Velocity)
- MV Vmean (MV Mean Velocity)
- MV VTI (Mitral Valve Velocity Time Integral)
- MVA (Mitral Valve Area)
- MVA By PHT (Mitral Valve Area according to PHT)

- MVA by plan (Mitral Valve Area, 2D)
- MVET (Mitral Valve Ejection Time)
- P Vein A (Pulmonary Vein Velocity Peak A) – reverse
- P Vein A Dur (Pulmonary Vein A-Wave Duration)
- P Vein D (Pulmonary Vein End-Diastolic Peak Velocity)
- P Vein S (Pulmonary Vein Systolic Peak Velocity)
- PAEDP (Pulmonary Artery Diastolic Pressure)
- PE(d) (Pericard Effusion, M-mode)
- PEs (Pericard Effusion, 2D)
- PR max PG (Pulmonic Insuf. Peak Pressure Gradient)
- PR mean PG (Pulmonic Insuf. Mean Pressure Gradient)
- PR PHT (Pulmonic Insuf. Pressure Half Time)
- PR Vmax (Pulmonic Insuf. Peak Velocity)
- PR VTI (Pulmonic Insuf. Velocity Time Integral)
- PRend max PG (Pulmonic Insuf. End-Diastole Pressure Gradient)
- PRend Vmax (Pulmonic Insuf. End-Diastolic Velocity)
- Pulmonic Diam (Pulmonary Artery Diameter, 2D)
- PV Acc Slope (Pulmonic Valve Flow Acceleration)
- PV Acc Time (Pulmonic Valve Acceleration Time)
- PV Acc Time/ET Ratio (PV Acceleration to Ejection Time Ratio)
- PV an diam (Pulmonic Valve Annulus Diameter, 2D)
- PV Ann Area (Pulmonic Valve Area)
- PV CO (Cardiac Output by Pulmonic Flow)
- PV max PG (Pulmonic Valve Peak Pressure Gradient)
- PV mean PG (Pulmonic Valve Mean Pressure Gradient)
- PV SV (Stroke Volume by Pulmonic Flow)
- PV Vmax (Pulmonary Artery Peak Velocity)
- PV Vmean (PV Mean Velocity)
- PV VTI (Pulmonic Valve Velocity Time Integral)
- PVA (VTI) (Pulmonary Artery Velocity Time Integral)
- PVein S/D Ratio (Pulmonary Vein SD Ratio)
- PVET (Pulmonic Valve Ejection Time)
- PVPEP (Pulmonic Valve Pre-Ejection Period)
- PVPEP/ET Ratio (PV Pre-Ejection to Ejection Time Ratio)
- Qp/Qs (Pulmonic-to-Systemic Flow Ratio)
- RA Major (Right Atrium Major, 2D)
- RA Minor (Right Atrium Minor, 2D)
- RAA (d) (Right Atrium Area, 2D, Diastole)
- RAA (s) (Right Atrium Area, 2D, Systole)
- RAEDV A2C (Right Atrium End Diastolic Volume, Apical 2 Chamber)
- RAESV A-L (RA End Systole Volume [A-L])
- RALd (Right Atrium Length, Diastole)
- RALs (RA Length, Systole)
- RIMP (Right Index of Myocardial Performance)
- RJA (A4C) (Regurgitant Jet Area)
- RJA/LAA (Regurgitant Jet Area ratio RJA/LAA)
- RV Major (Right Ventricle Major)
- RV Minor (Right Ventricle Minor)
- RV S' (Tricuspid annulus systolic excursion velocity)
- RVAWd (Right Ventricle Wall Thickness, Diastolic, 2D)
- RVAWs (Right Ventricle Wall Thickness, Systolic, 2D)
- RVET (Right Ventricle Ejection Time)
- RVIDd (Right Ventricle Diameter, Diastolic, 2D)
- RVIDs (Right Ventricle Diameter, Systolic, 2D)
- RVOT Area (Right Ventricle Outflow Tract Area)
- RVOT Diam (RV Output Tract Diameter, 2D)
- RVOT Diam (RV Output Tract Diameter, M-Mode)
- RVOT max PG (RVOT Peak Pressure Gradient)
- RVOT meanPG (RVOT Mean Pressure Gradient)
- RVOT SI (LV Stroke Volume Index by Pulmonic Flow)
- RVOT SV (Stroke Volume by Pulmonic Flow)
- RVOT Vmax (RVOT Peak Velocity)
- RVOT Vmean (RVOT Mean Velocity)
- RVOT VTI (RVOT Velocity Time Integral)
- RVSP (Right Ventricle Systolic Pressure)
- RVWd (Right Ventricle Wall Thickness, Diastolic, M-mode)
- RVWs (Right Ventricle Wall Thickness, Systolic, M-mode)
- RAA (d) (Right Atrium Area, 2D, Diastole)
- RAA (s) (Right Atrium Area, 2D, Systole)
- SI (A-L A2C) (LV Stroke Index, Single Plane, 2CH, Area-Length)
- SI (A-L A4C) (LV Stroke Index, Single Plane, 4CH, Area-Length)
- SI (Bi-plane) (LV Stroke Index, Bi-plane, MOD)
- SI (bullet) (LV Stroke Index, Bi-plane, Bullet)
- SI (MOD A2C) (LV Stroke Index, Single Plane, 2CH, MOD)
- SI (MOD A4C) (LV Stroke Index, Single Plane, 4CH, MOD)
- SI (Teich) (LV Stroke Index, Teicholtz, 2D)
- SI (Teich) (LV Stroke Index, Teicholtz, M-mode)
- SV (A-L A2C) (LV Stroke Volume, Single Plane, 2CH, Area-Length)
- SV (A-L A4C) (LV Stroke Volume, Single Plane, 4CH, Area-Length)
- SV (Bi-plane) (LV Stroke Volume, Bi-plane, MOD)

- SV (bullet) (LV Stroke Volume, Bi-plane, Bullet)
- SV (MOD A2C) (LV Stroke Volume, Single-plane, 2CH, MOD) – Simpson
- SV (MOD A4C) (LV Stroke Volume, Single-plane, 4CH, MOD) – Simpson
- SV (Cube) (LV Stroke Volume, 2D, Cubic)
- SV (Cube) (LV Stroke Volume, M-mode, Cubic)
- SV (Teich) (LV Stroke Volume, 2D, Teicholtz)
- SV (Teich) (LV Stroke Volume, M-mode, Teicholtz)
- Systemic Diam (Systemic Vein Diameter, 2D)
- Systemic Vmax (Systemic Vein Peak Velocity)
- Systemic VTI (Systemic Vein Velocity Time Integral)
- TAPSE (Tricuspid Annular Plane Systolic Excursion)
- TCO (Tricuspid Valve Closure to Opening)
- TR max PG (Tricuspid Regurg. Peak Pressure Gradient)
- TR mean PG (Tricuspid Regurg. Mean Pressure Gradient)
- TR Vmax (Tricuspid Regurg. Peak Velocity)
- TR Vmean (Tricuspid Regurg. Mean Velocity)
- TR VTI (Tricuspid Regurgitation Velocity Time Integral)
- TV A dur (Tricuspid Valve A-Wave Duration)
- TV A Velocity (Tricuspid Valve A Velocity)
- TV Acc Time (Tricuspid Valve Time to Peak)
- TV Ann Area (Tricuspid Valve Area)
- TV ann diam (Tricuspid Valve Annulus Diameter, 2D)
- TV Area (Tricuspid Valve Area, 2D)
- TV CO (Cardiac Output by Tricuspid Flow)
- TV Dec Slope (Tricuspid Valve Flow Deceleration)
- TV E Velocity (Tricuspid Valve E Velocity)

- TV E/A Ratio (Tricuspid Valve E-Peak to A-Peak Ratio)
- TV max PG (Tricuspid Valve Peak Pressure Gradient)
- TV mean PG (Tricuspid Valve Mean Pressure Gradient)
- TV mean PG (Tricuspid Valve Mean Pressure Gradient)
- TV PHT (Tricuspid Valve Pressure Half Time)
- TV SV (Stroke Volume by Tricuspid Flow)
- TV Vmean (TV Mean Velocity)
- TV VTI (Tricuspid Valve Velocity Time Integral)
- VSD max PG (VSD Peak Pressure Gradient)
- VSD Vmax (VSD Peak Velocity)

Please refer to the Reference Manual for the full list of measurements and calculations for all applications.

Annotations

Body Marks

- Body mark icons for location and position of probe
- Option to automatically activate body mark on freeze
- Easy selection of body marks from touch screen
- Easy selection of body marks for dual-screen layout

Text Annotations

- Easy selection of text annotations from touch screen
- Option to automatically activate annotation on freeze

Scan Assist Pro (optional)

- Customizable automations that assist the user through each step of the scan
- Facilitates consistency and reduce keystrokes
- Supports selection of all modes, all measurements and dual annotations

- Imaging attributes: Octave, Steer, Dual/Quad screen, Compound, LOGIQ View, Zoom, Depth, Scale and Baseline
- On-line or off-line protocol editor
- Image acquisition according to predefined protocol templates
- Various factory protocol templates
- User-configurable protocol templates

Stress Echo (optional)

Supported Protocol Examinations

- 2D pharmacological stress echo
- 2D bicycle stress echo
- 2D continuous capture stress echo (treadmill stress echo)
- Cardiac resynchronization therapy programming protocols (available with the Advanced QScan option)

Protocol Examinations Features (enabled with stress option)

- Wall motion scoring: Analysis by wall motion in individual myocardial segments
- Show reference: Show a reference image from baseline or previous level during acquisition
- Smart stress: Automatically set up various scanning parameters (for instance geometry, frequency, gain, etc.) according to same projection on previous level
- Scan mode settings: Scan mode may be specified for individual views in the protocol
- Preview of store: Show running loops as preview before storing to the examination

Continuous Capture

- Continuously acquire large amounts of 2D image data, and selection of projection views for analysis afterwards
- The entire continuous capture recording may be kept in memory while it is possible to store new images outside the protocol template, or the entire recording can be stored to file
- Selection of projection views on scanner or EchoPAC when the entire recording is stored to file

Wall Motion Scoring

- As part of the measurement and analysis package one can access a wall motion assessment module, providing analysis/scoring of individual myocardial segments
- For use with all stress modalities

Cardiac Resynchronization Therapy (CRT) Programming Protocols

- CRT protocols require Stress and Advanced QScan
- Tailored acquisition protocol for data needed for programming of AV and VV delays in biventricular pacemakers
- Image acquisition of a set of projection views with various scan mode settings
- Template editor
- User-configurable protocol templates
- Configure protocol name, number of levels and views, name of level and views and several other protocol settings (smart stress, show reference, scan mode, preview of store, timer handling, etc.)

CARTO[®] 3 Interface (optional)

- The system can interface with the CARTO 3 EP navigation system and the SOUNDSTAR[®] ultrasound catheters manufactured by Biosense Webster, Inc
- The interface will allow the Vivid S60N system to send images to the CARTO 3 EP system over a video cable
- The Vivid S60N is able to send ultrasound scaling parameters to the CARTO 3 EP system via a peer-to-peer LAN connection

Safety Conformance

- The Vivid S60N is built to meet the requirements of:
 - IEC60601-2-37
 - IEC60601-1
 - IEC60601-1-2
 - IEC60601-1-6
 - UL60601-1
 - NEMA UD3

- The European Medical Devices Directive (MDD) 93/42/EEC (CE Mark)
- Directive 2011/65/EU on the restriction of use of certain hazardous substances
- The Vivid S60N ultrasound unit is a Class I device, with BF (probes) and CF (ECG leads) applied parts according to IEC60601-1
- The Vivid S60N ultrasound unit meets the EMC requirements in IEC/EN60601-1-2:2007 Class B

Security

Virus Protection

To reduce virus vulnerability, Vivid S60N is configured with a minimal set of open ports and with all network services not actively used by the system closed down. This helps to significantly reduce the risk of a virus attack on Vivid S60N.

GE is continuously judging the need for additional actions to reduce vulnerability of equipment; this includes vulnerability scanning of our products and evaluation of new security patches for the 3rd-party technology used. Microsoft[®] (and other) security patches that address serious issues with Vivid S60N will be made available to customers after GE verification of those patches.

Whitelisting

- Prevents non-listed applications from running

User Policies

- Secure and advanced user password and login scheme according to user's password requirements

LDAP

- Users can log in to the system by using the same user credentials as used for domain connected computers

Disc Encryption

- Optional encryption of the scanner's E drive containing patient identifiable data

Transducers

3Sc-RS Phased Array Probe

- Probe presets: Cardiac adult, pediatric, abdomen, fetal heart, adult transcranial, coronary, stress, LVO contrast, LVO stress, OB/GYN, vascular
- Biopsy guide: Multi-angle disposable with a reusable bracket

6S-D Phased Array Probe

- Probe presets: Pediatric, cardiac, coronary, neonatal head, fetal heart, abdominal

12S-D Phased Array Probe

- Probe presets: Pediatric, coronary, neonatal head, vascular (incl. carotid, LEA, LEV, UEA, UEV), abdomen

9L-D Linear Array Probe

- Probe presets: Vascular (incl. carotid, LEA, LEV, UEA, UEV), musculoskeletal conventional, nerves, small parts, thyroid
- Biopsy guide: Multi-angle disposable with a reusable bracket

11L-D Linear Array Probe

- Probe presets: Vascular, small parts, breast, thyroid, scrotal, musculoskeletal conventional, nerves
- Biopsy guide: Multi-angle disposable with reusable bracket

L8-18i-D Intraoperative Linear Array Probe

- Probe presets: Cardiac, musculoskeletal, vascular, small parts

C1-5-D Curved Array Probe (Convex)

- Probe presets: Abdomen, OB/GYN, urology, vascular, fetal heart
- Biopsy guide: Multi-angle, disposable with a reusable bracket

C1-6-D XDclear Single Crystal Curved Array Probe (Convex)

- Probe presets: Abdomen, OB/GYN, urology, vascular, fetal heart
- Biopsy guide: Multi-angle, disposable with a reusable bracket

C2-9-D XDclear

Single Crystal Curved Array Probe (Convex)

- Probe presets: Abdomen, OB/GYN, urology, fetal heart
- Biopsy guide: Multi-angle disposable with reusable bracket

C3-10-D XDclear

Single Crystal Curved Array Probe (Micro Convex)

- Probe presets: Neonatal head, vascular, abdomen, musculoskeletal, nerves

IC5-9-D Convex (Endocavity) Probe

- Probe presets: OB/GYN, urology, fetal heart
- Biopsy guide: Single angle, disposable bracket

P2D Pencil Probe

- Probe preset: Cardiac

P6D Pencil Probe

- Probe preset: Vascular

6Tc-RS TEE Probe

- Probe presets: Cardiac, coronary, LVO contrast

6VT-D Active Matrix TEE Probe

- Probe presets: Cardiac, coronary, LVO contrast

9T-RS TEE Probe

- Probe preset: Pediatric cardiac

Catheter Cable ICE Probe Connector

- Allows connecting the AcuNav® ICE catheters to Vivid S60N

ACUSON® AcuNav 10F IntraCardiac Echo (ICE) Catheter⁴

- Probe presets: ICE

ACUSON AcuNav 8F IntraCardiac Echo (ICE) Catheter⁴

- Probe presets: ICE

SOUNDSTAR 3D Ultrasound Catheter based on AcuNav 10F IntraCardiac Echo (ICE) Catheter⁴

- Probe presets: ICE, CARTO

SOUNDSTAR eco 10F G Ultrasound Catheter based on AcuNav 10F IntraCardiac Echo (ICE) Catheter⁴

- Probe presets: ICE, CARTO

SOUNDSTAR eco 8F G Ultrasound Catheter based on AcuNav 8F IntraCardiac Echo (ICE) Catheter⁴

- Probe presets: ICE, CARTO

Wideband Probes

- Electronic selection between four solid-state and one stand-alone Doppler probe connectors
- Three probe sockets are DLP type
- One RS probe socket for support of TEE and ICE probes

PROBE	FREQUENCY RANGE	CATALOG #
3Sc-RS (Sector)	1.3 – 4.5 MHz	H45041DL
6S-D (Sector)	2.4 – 8.0 MHz	H45021RR
12S-D (Sector)	4.0 – 12.0 MHz	H45021RT
9L-D (Linear)	2.4 – 10.0 MHz	H40442LM
11L-D (Linear)	4.0 – 12.0 MHz	H40432LN
L8-18i-D (Intraoperative Linear)	5.0 – 18.0 MHz	H40452LL
C1-5-D (Convex)	1.4 – 6.0 MHz	H40452LE
C1-6-D (Convex)	1.4 – 6.0 MHz	H40472LT
C2-9-D (Convex)	2.3 – 8.4 MHz	H40462LN
C3-10-D (Micro Convex)	3.0 – 10.0 MHz	H40482LB
IC5-9-D (Convex Endocavity)	3.3 – 8.6 MHz	H40442LK
P2D (Pencil)	2.0 MHz	H4830JE
P6D (Pencil)	6.3 MHz	H4830JG
6Tc-RS (TEE)	3.0 – 8.0 MHz	H45551ZE
6VT-D (TEE) ⁵	3.0 – 8.0 MHz	H45581BJ
9T-RS (TEE)	3.0 – 10.0 MHz	H45531YM
Catheter Cable ICE probe connector		H48952AR
ACUSON AcuNav 10F ⁴	4.5 – 11.5 MHz	Distributed by Biosense Webster, Inc
ACUSON AcuNav 8F ⁴	4.5 – 11.5 MHz	Distributed by Biosense Webster, Inc
SOUNDSTAR 3D Ultrasound Catheter based on AcuNav 10F ⁴	4.5 – 11.5 MHz	Distributed by Biosense Webster, Inc
SOUNDSTAR eco 10F G Ultrasound Catheter based on AcuNav 10F ⁴	4.5 – 11.5 MHz	Distributed by Biosense Webster, Inc
SOUNDSTAR eco 8F G Ultrasound Catheter based on AcuNav 8F ⁴	4.5 – 11.5 MHz	Distributed by Biosense Webster, Inc

⁴ Not available in all countries. Please contact Biosense Webster, Inc for availability.

⁵ Also 6VT-D with catalog # H45561TA is supported.

Product may not be available in all countries and regions.
Full product technical specification is available upon request.
Contact a GE Healthcare Representative for more information.
Please visit www.gehealthcare.com/promotional-locations.

Data subject to change.

© 2018 General Electric Company. DOC2128811

GE, the GE Monogram, imagination at work, Vivid, XDclear, EchoPAC, LOGIQ, InSite and cSound are registered trademarks of General Electric Company.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

Windows and Microsoft are registered trademarks of Microsoft Corporation.

Sony is the registered trademark of Sony Corporation.

Tricefy is a trademark of Trice Imaging, Inc.

ACUSON and AcuNav are trademarks of Siemens Medical Solutions USA, Inc.

CARTO and SOUNDSTAR are trademarks of Biosense Webster, Inc.

All other third-party trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

About GE Healthcare

GE Healthcare is a leading provider of medical imaging, monitoring, biomanufacturing, and cell and gene therapy technologies. GE Healthcare enables precision health in diagnostics, therapeutics and monitoring through intelligent devices, data analytics, applications and services. With over 100 years of experience and leadership in the healthcare industry and more than 50,000 employees globally, GE Healthcare helps healthcare providers, researchers and life sciences companies in their mission to improve outcomes for patients around the world. Follow us on [Facebook.com](https://www.facebook.com/gehealthcare), [LinkedIn.com](https://www.linkedin.com/company/gehealthcare), [Twitter.com](https://twitter.com/gehealthcare) and [The Pulse.com](http://ThePulse.com) for latest news, or visit our website www.gehealthcare.com for more information.

GE Healthcare
9900 Innovation Drive
Wauwatosa, WI 53226
U.S.A.

www.gehealthcare.com



imagination at work