

EVO

NEW

THE



TAKE ANOTHER LOOK...

For nearly 40 years, we've manufactured animal-specific ultrasound right here in the US. Now we've done it again... the next generation of our flagship system, the **NEW EVO**.

With new high-resolution probes, breakthrough imaging and battery performance, the **EVO III** provides high-quality ultrasound in a rugged, portable machine.

- New broadband probes for all species
- Completely redesigned image software with the latest technology
- Customized, veterinary-specific exam presets
- Sunlight-readable LED display
- WiFi streaming via EVOStream app
- Ruggedized case, splash resistant, and portable [just over 6 lbs]—sealed keyboard and probe connections allow easy cleaning
- Interchangeable batteries, for continuous use
- Designed and manufactured in Colorado

E.I. Medical Imaging[®]
1.866.365.6596 • www.eimedical.com/the_new_evo

TRANSDUCERS NEW, IMPROVED!!



eCL3
2-5 MHz convex transducer. Abdominal, reproduction, lung in large animals and exotics.



eCLi4
3-7 MHz broadband array. Curved rectal probe for bovine, equine reproduction; abdomen, exotics, small ruminants.



eL14
6-14 MHz broadband array for equine MSK, small animal, exotics.




eC6
4-9 MHz broadband array, 18 cm penetration. Small animal abdominal, exotics.



eC9
5-10MHz broadband array. Cats, small dogs, exotics, OPU.



eL7
5-9 MHz broadband rectal linear for reproductive exams. Bovine and equine reproduction, tendon and lung.



L3ASE
Linear transducer, Bovine carcass merit.



eC90PU
5-10 MHz broadband array for bovine and equine ovum pickup. Adjustable handle, small height, lightweight.

PROBE ACCESSORY



I.C.E.®
The IBEX Customizable Extension lessens fatigue and stress on the user's arm and shoulder. New, improved version is 25% lighter, smaller diameter, more durable than previous models. The I.C.E. extension allows easy insertion or removal of the probe without damaging cables. Compatible with all IBEX and EVO rectal probes.

A WORD ABOUT OUR CABLES...

IBEX® cables tested to 1.2 million flex cycles without breaking

- Custom designed cables are abrasion resistant
- Only manufacturer to design all transducers specifically for animal use
- Kevlar reinforced

...AND OUR NEW TRANSDUCERS

As part of our commitment to innovation, we have re-designed our probes for improved penetration and resolution. Broadband, variable frequency, our transducers are designed specifically for rugged environments. But if the worst happens, they can be repaired right here in our Colorado facility.

ACCESSORIES

i3 LED goggles, EVO sun shade, Bluetooth remote, EVO waist pouch—all designed to streamline your workflow! www.eimedical.com/accessories

E.I. MEDICAL IMAGING

E.I. Medical Imaging® was founded in 1984 with the vision of developing state-of-the-art, real-time ultrasound to serve veterinarians and livestock producers world-wide.

Our products have evolved with the needs of the market. We are the only manufacturer of portable ultrasound systems engineered and developed in the United States specifically for the animal industry.



CONTACT US

1.866.365.6596
970.669.1793
970.669.1902 [fax]
info@eimedical.com