

TRANSDUCERS/PROBES

EVO



eCL3
2-5 MHz convex transducer. Abdominal, reproduction, lung in large animals and exotics.



eCLi4
3-7 MHz broadband array, 25 cm depth. Curved rectal probe for bovine, equine reproduction; abdomen, exotics, small ruminants.



eL7
5-9 MHz broadband rectal linear for reproductive exams. Bovine and equine reproduction, tendon and lung.



eC6
4-9 MHz broadband array, 18 cm penetration. Small animal abdominal, exotics.



eL14
6-14 MHz broadband array for equine MSK, small animal, exotics.



eC9
5-10 MHz broadband array. Cats, small dogs, exotics, OPU.



L3ASE
Linear transducer, Bovine carcass merit.



eC90PU
5-10 MHz broadband array for bovine and equine ovum pickup. Adjustable handle, small height, lightweight.

PRO, LITE, LITE HR



iL7
5-9 MHz broadband rectal linear for reproductive exams. Bovine and equine reproduction, tendon and lung.



iCL3
Low-frequency [2.5-5 MHz] 24 cm depth, deep penetrating convex transducer. Large animal abdomen, exotics, small ruminant.



iCLi4
Side fire curved [2.5-5 MHz], 24 cm depth. Bovine, equine reproduction; large animal abdomen, exotics, small ruminants.



L3.8AS
Linear transducer, Bovine carcass merit.



iC90PU
5-10 MHz broadband array for bovine and equine ovum pickup. Adjustable handle, small height, lightweight.

PROBE ACCESSORY



I.C.E.® NEW AND IMPROVED!
The IBEX Customizable Extension lessens fatigue and stress on the user's arm and shoulder. New, improved version is 25% lighter, smaller diameter, more durable than previous models. The I.C.E. extension allows easy insertion or removal of the probe without damaging cables. Compatible with all IBEX and EVO rectal probes.

A word about our cables...

IBEX® cables tested to 1.2 million flex cycles without breaking

- Custom designed cables are abrasion resistant
- Only manufacturer to design all transducers specifically for animal use
- Kevlar reinforced

...and our transducers

All of our transducers have been re-designed and upgraded!

As part of our commitment to innovation, we have re-designed our probes for improved penetration and resolution.

Broadband, variable frequency, our transducers are designed specifically for rugged environments. But if the worst happens, they can be repaired right here in our Colorado facility.

TRANSDUCERS/PROBES

SA



eCL3
2-5 MHz convex transducer. Abdominal, repro, lung in large animals and exotics.



eC9
5-10MHz broadband array. Cats, small dogs, exotics, OPU.



eC6
4-9 MHz broadband array, 18 cm penetration. Small animal abdominal, exotics.



eL14
6-14 MHz broadband array for equine MSK, small animal, exotics.

i3 LED GOGGLES

The i3 LED goggles are small, lightweight and compatible with all IBEX® veterinary ultrasound models as well as with Bantam and Easi-Scan [models 1-4].

They come with versatile mounting options: ratcheted headband, head strap and hat clip. Also included—detachable glare shields to block sunlight.

Custom designed optics

- LCD display with LED backlight
- Half the weight of InSite2
- Specially designed eye cups for inside and outdoor viewing
- GoPro mount
- Standard, extended warranties available



ACCESSORIES



WiFi Video Transmitter

- Connects to all IBEX PRO and LITE ultrasound series
- Streams ultrasound image to app on tablet or phone
- Creates private network with no WiFi service required
- Uses goggles port—splitter is available if goggles will be used simultaneously



IBEX Waist Pouch

- Accommodates the IBEX LITE and PRO series [in closed position]
- Keeps unit self-contained and streamlined when scanning with headset



EVO Remote

- Bluetooth connection—auto connects once paired with an EVO ultrasound
- 7 keys can be reassigned to any key for versatility
- Use to control gain, depth, save and record images, or other custom features



EVO Waist Pouch

- Accommodates EVO machines [in closed position]
- Keeps unit self-contained and streamlined when scanning with headset



EVO Sun Shade

- Magnetic shade snaps on to EVO screens to allow easy viewing in bright conditions

CONTACT US

E.I. Medical Imaging
815 14th St. SW, Unit C210
Loveland, CO 80537
www.eimedical.com/get-info

toll-free/1.866.365.6596
phone/970.669.1793
fax/970.669.1902
info@eimedical.com

