

IBEX[®] LITE

portable | rugged | ultrasound



E.I. Medical Imaging[®]



FEATURES

- DuraScan® case: water, dust and impact resistant
- 3.3 lb [1.5 kg]
- B, B/M modes
- Auto enhanced field of view [auto adjusts focus, frequency, gain and other parameters]
- Digital beam-former
- Gestational tables
- Stores still frames and video loops
- Li-Ion smart batteries with up to 6 hours of run time
- Multiple headset options available
- I.C.E. arms-free scanning available
- Linear, convex, micro-convex and specialty transducers available
- Standard, extended warranties available

LITE

With 90% of the features and the same image quality of the IBEX® PRO, the IBEX LITE is the ideal choice for bovine ultrasound. From basic repro to fetal sexing and calf lungs, the IBEX LITE has all the features you need in an ultra light, ruggedized package.

Without the added weight and cost of the screen, the IBEX LITE is ideal for dairy practice and outdoor field work. The video headset is great for working in both bright sunlight and low light conditions, and there are a variety of interchangeable transducers available. Worn around the waist, the unit is lightweight and makes for streamlined workflow. Like the other products in the IBEX line, the LITE is rugged and water-resistant, with long battery life.



TRANSDUCERS [LITE ONLY]



L8.0
High-frequency linear transducer
Equine tendons, small animal abdomens, neurovascular applications, superficial structures



MC8.0
High-frequency microconvex transducer
Small animal abdomen, bovine reproduction, small ruminant pregnancies... used in the OPU aspiration device



CL3Ei
Low-frequency, deep penetrating convex transducer
Large animal abdomen, exotics, small ruminant



I.C.E.® NEW AND IMPROVED
The *updated* IBEX® Customizable Extension lessens fatigue and stress on user's arm and shoulder—25% lighter, smaller diameter, more durable than previous models
The I.C.E. extension allows easy insertion or removal of the probe without damaging cables
Compatible with all IBEX and EVO rectal probes



CL4HDi
Side fire curved
Bovine reproduction, large animal abdomen, exotics, small ruminant



L7HDi
Linear transducer
Bovine and equine reproduction, tendon



L3.8AS
Linear transducer
Bovine carcass merit



C90PUi
Ovum pickup probe and handle [5-10 MHz], 14 cm depth
Adjustable handle, shallower focus for better imaging; small height, lightweight
Improved positioning of transducer array and needle exit
Equine and bovine

ACCESSORIES



i3 LED GOGGLES
The *i3 LED goggles* are small, lightweight and compatible with all IBEX veterinary ultrasound models
Versatile mounting options: ratcheted headband, head strap and hat clip.
Comes with detachable glare shields to block sunlight.
LCD display with LED backlight with specially designed eye cups for inside/outdoor viewing
Half the weight of InSite2
GoPro mount
Standard, extended warranties available



IBEX Waist Pouch
Accommodates the IBEX LITE and PRO series [in closed position]
Keeps unit self-contained and streamlined when scanning with headset



WiFi Video Transmitter
Connects to all IBEX PRO and LITE veterinary ultrasound series
Streams ultrasound image to app on tablet or phone
Creates private network with no WiFi service required
Uses goggles port—splitter is available if goggles will be used simultaneously

Transducers + cables

- All IBEX cables tested to 1.2 million flex cycles:
 - Custom designed cable jacket is extremely durable, chemical resistant and bio-compatible
 - Designed specifically for use in the tough animal environment
 - Kevlar strength members embedded in cable bundle to protect coax breakage
 - All materials RoHS compliant
- All transducers are designed specifically for rugged environments. But if the worst happens, they can be repaired, too!

E.I. MEDICAL IMAGING®

E.I. Medical Imaging was founded in 1984 with the vision of developing a state-of-the-art, real-time ultrasound scanner to serve veterinarians and livestock producers world-wide.

Our products have evolved with the needs of the market. We are proud to be the only manufacturer of portable ultrasound systems engineered and developed in the United States specifically for the animal industry.

For over 30 years, the company's core values have remained intact: putting the customer first and delivering solid, effective ultrasound solutions. Today we are a world leader in the durable, portable, veterinary ultrasound industry.

CONTACT US

E.I. Medical Imaging

110 12th St. SW, Unit 102

Loveland, CO 80537-6396

www.eimedical.com/get-info

toll-free/1.866.365.6596

phone/970.669.1793

fax/970.669.1902

info@eimedical.com



EDUCATION

EIMI is committed to education. From online learning resources to classroom and wet lab training experiences, we are here to ensure that you get the most out of your ultrasound equipment. Our on-staff veterinarian and sales team are always available to assist you as you put your IBEX® system to work.

EIMI is also involved in educating future veterinarians, and offers two specific programs toward that end:

- EIU™ Student Externship Program—EIMI provides ultrasounds on loan to vet students for externship experiences. Contact us for more information.
- EIU™ Seminars—EIMI hosts ultrasound seminars at veterinary schools and in the field, comprising a short introductory lecture and hands-on wet labs. Contact us to schedule.

CUSTOMER SUPPORT/LOANER PROGRAM

Downtime is not in the vernacular at E.I. Medical Imaging®... service is the hallmark of our business. We engineer systems specifically for veterinary use and understand the challenges downtime brings. With industry-leading service turnaround times and a loaner program available, we ensure continuity of your business. All service is conducted at our EIMI facility in Loveland, Colorado, using only genuine EIMI parts.