GE HealthCare

Versana Ultrasound

Veterinary care with confidence





Versana Essential™



Versana Premier™



Versana Balance™



Versana Active™

Images that speak for your patient

We designed our Versana ultrasound systems to provide high-quality images that can help you provide better veterinary care. Backed by GE HealthCare's legacy of ultrasound excellence, Versana ultrasound systems provide world-class imaging and offers broad applicability to support comprehensive care for a wide variety of patients.

No matter which Versasna system you choose, you'll be able to conduct a diverse array of exams, plus you'll have access to clinical tools that expedite workflows and simplify reporting. And, if you're new to ultrasound, rest assured that Versana ultrasound systems are very easy to learn. These systems are easy to use and backed with robust GE HealthCare service and support.

- **World-class image quality**
- **Easy to learn**
- **Broad applicability**
- **Backed with service and support**

A complete solution

Improve care with broad clinical capabilities

We know the cases veterinarians see daily can vary greatly, so Versana ultrasound systems are designed to help handle whatever comes your way. Our systems feature an intuitive design, user-friendly interface, and simple workflows so you easily supplement exams with ultrasound imaging whenever you feel the need. Versana ultrasound systems are so versatile, you can use them to perform many types of exams such as:

- Cardiac
- Abdominal
- Pelvic
- Musculoskeletal
- Thyroid
- Obstetrics
- Vascular
- Ultrasound-guided procedures



Clinical applications[1] that simplify exams :

3D/4D Capability
Aids decision making by enabling additional dimensions in real time

Needle Recognition
Simplifies biopsy procedures with technology that clarifies the exact position of the needlepoint

Whizz B-mode
Allows you to continuously optimize brightness, contrast, gain and uniformity of images while scanning different tissues

Whizz color
Optimizes by adjusting Color Flow gain and frequency automatically when activated by user

Whizz Pulsed/Continuous wave Doppler
Automatically adjusts the velocity scale/PRF (live imaging only), baseline shift, and invert (if preset)

Elastography
Shows the spatial distribution of tissue elasticity properties in a region of interest by estimating the strain before and after tissue distortion caused by external and internal forces

Follow-up Tool
Enables comparison of previous and current exams simultaneously

A complete solution

Designed to help expedite exams

Although ultrasound is painless, keeping animals in the right position can require some effort. For this reason, Versana ultrasound systems are designed with tools[1] you can use to perform scans as quickly as possible.

Scan Assistant
Enables exam protocols that walk users through all steps of image acquisition, saving keystrokes, reducing stress and fatigue, and maximizing exam consistency

Patient information tool
Shortens exam set-up time by including fields to log animal ID, animal name, host, species, and weight


OB measurement options
Speeds exams with editable and customizable fields

Animal specific body markers and annotations
Creates efficiencies with animal specific body markers and annotations


Animal-specific report templates
Speeds reporting with OB templates for cat, dog, equine; and create your own with template designer




Microconvex imaging Liver



Sector Cardiac imaging Color Doppler



Linear with Virtual Convex imaging Spleen



Linear imaging Renal

You'll also appreciate these additional features that help simplify system maintenance and accelerate exams:

- Fast boot-up
- Large frame high resolution LED backlit monitor display[2]
- Veterinary probes
- Cable management
- GE HealthCare's latest cybersecurity updates
- Large internal storage[2]
- Probe check to ensure optimal probe performance



Opportunity for growth

Meet your needs now and in the future

One of the best things about Versana ultrasound systems is that they are designed to meet your needs today and grow with you tomorrow. GE HealthCare is there to support you during and after your purchase with in-depth training and education.[3] Additionally, My Trainer, a service training module on the system, will help you get started quickly.



Flexibility to adapt as you grow

As your clinical confidence grows, you can expand capabilities with a range of transducers and applications that enhance the amount of information you see during patient exams.

And finally, you can always upgrade your current Versana Ultrasound system or upgrade to a different GE HealthCare ultrasound system as your practice advances.

Partner you can trust

Quality and reliability guaranteed

Ultrasound can quickly become a highly-used tool in your practice, and that's just one more reason to trust Versana ultrasound systems. Years of research have gone into the development of our ultrasound technology, and you'll find our best-in-class image quality is a perfect way to take your practice to the next level.

Easy to update and keep secure

You don't need to be an IT expert to use Versana ultrasound systems. Features like e-delivery[4] of software updates help keep systems current and advanced cybersecurity and data protections ease the burden of guarding patient information.

Rugged and dependable

Because animal care equipment can take a beating, Versana ultrasound systems undergo rigorous testing to ensure long performance. Our systems air filter prevents fur from being drawn into the central processing unit (CPU). In the event you require service or support, we offer remote and on-site options plus service engineers are always close by should you need them. For added peace-of-mind, flexible, multi-year service coverage and warranties are available to help you deliver quality, reliable care, year after year. And of course, we offer flexible financing options[5] which means there's never been a better time to bring a Versana ultrasound system onboard. Check with your GE HealthCare provider to learn which Versana ultrasound system is right you.



Versana ultrasound systems are virtually worry-free:

- E-delivery[4] of software updates
- Advanced cybersecurity
- Remote and on-site support options
- Multi-year service coverage and warranties
- Flexible financing options

GE HealthCare

About GE HealthCare

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 100 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from prevention and screening, to diagnosis, treatment, therapy, and monitoring. We are an $18 billion business with 51,000 employees working to create a world where healthcare has no limits.

Follow us on Facebook, LinkedIn, Twitter, Instagram and Insights for the latest news, or visit our website gehealthcare.com for more information.

References:

1. Clinical functions may not be available on all probes and are revision dependent. Please check with your local GE HealthCare representative for more details and availability.
2. Size and/or capacity varies depending on the Versana Ultrasound System.
3. Local education offerings may vary. Please check with your local GE HealthCare representative.
4. Available on Versana Balance and Versana Essential. Not available on Versana Active and Versana Premier.
5. Local educational options may vary according to region. Please check with your local GE HealthCare representative.

© 2023 GE HealthCare – All rights reserved.

GE HealthCare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE HealthCare representative for the most current information. Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Versana Active, Versana Balance, Versana Essential, and Versana Premier are trademarks of GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

JB25088XX